
Chapter 3- Precision Cutting

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WAFERING BLADE SELECTION GUIDELINES

Material	Characteristic	Speed (rpm)	Load (grams)	Blade (grit/conc.)
Silicon substrate	soft/brittle	<300	<100	Fine/low
Gallium arsenide	soft/brittle	<200	<100	Fine/low
Boron composites	very brittle	500	250	Fine/low
Ceramic fiber composites	very brittle	1000	500	Fine/low
Glasses	brittle	1000	500	Fine/low

Minerals	friable/brittle	>1500	>500	Fine/low
Alumina ceramic	hard/tough	>1500	>500	Medium / low
Zirconia (PSZ)	hard/tough	>3500	>800	Medium/low
Silicon nitride	hard/tough	>3500	>800	Medium/low
Metal matrix composites		>3500	>500	Medium/high
General purpose		variable	Variable	Medium/high

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Wafering Blades Description	Characteristic
Fine grit	10-20 micron (600 grit)
Medium grit	60-70 micron (220 grit)
Coarse grit	120 micron (120 grit)
High concentration	100%
Low concentration	50%

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

(a) Electroplated diamond wafering blade for cutting soft materials (left) and (b) alumina resin-rubber blade (right)

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




Partially stabilized zirconia sectioned at low speeds and low loads

Partially stabilized zirconia sectioned at high speeds and high loads

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3" Diameter (76.2 mm) Wafering Blade, 1/2" (12.7 mm) Arbor

Application	Conc.	Diamond size	Thickness	Catalog Number	Figure
Brittle materials (microelectronic materials, friable ceramics, minerals)	Low	Fine	0.006"	WB-0030LC	
Hard/tougher ceramics (Si ₃ N ₄ , ZrO ₂)	Low	Medium	0.006"	WB-0035LC	
Most metal samples	High	Medium	0.006"	WB-0035HC	

4" Diameter (102 mm) Wafering Blade, 1/2" (12.7 mm) Arbor

Application	Conc.	Diamond size	Thickness	Catalog Number	Figure
Brittle materials (microelectronic materials, friable ceramics, minerals)	Low	Fine	0.012"	WB-0040LC	
Hard/tougher ceramics (Si ₃ N ₄ , ZrO ₂)	Low	Medium	0.012"	WB-0045LC	
Most metal	High	Medium	0.012"	WB-0045HC	

samples					
---------	--	--	--	--	--

5" Diameter (127 mm) Wafering Blade, 1/2" (12.7 mm) Arbor

Application	Conc.	Diamond size	Thickness	Catalog Number	Figure
Brittle materials (microelectronic materials, friable ceramics, minerals)	Low	Fine	0.015"	WB-0050LC	
Hard/tougher ceramics (Si ₃ N ₄ , ZrO ₂)	Low	Medium	0.015"	WB-0055LC	
Most metal samples	High	Medium	0.015"	WB-0055HC	

6" Diameter (153 mm) Wafering Blade, 1/2" (12.7 mm) Arbor

Application	Conc.	Diamond size	Thickness	Catalog Number	Figure
Brittle materials (microelectronic materials, friable ceramics, minerals)	Low	Fine	0.020"	WB-0060LC	
Hard/tougher ceramics (Si ₃ N ₄ , ZrO ₂)	Low	Medium	0.020"	WB-0065LC	
Most metal samples	High	Medium	0.020"	WB-0065HC	

Electroplated Diamond Wafering Blades, 1/2" (12.7 mm) Arbor

Application	Blade Diameter	Catalog Number	Figure
-------------	----------------	----------------	--------

4-inch Electroplated Diamond blade for soft gummy materials such as plastics, rubber, PCB's and other soft materials	4-inch	EPD-0040	
5-inch Electroplated Diamond blade for soft gummy materials such as plastics, rubber, PCB's and other soft materials	5-inch	EPD-0040	
6-inch Electroplated Diamond blade for soft gummy materials such as plastics, rubber, PCB's and other soft materials	6-inch	EPD-0040	


CBN Wafering Blades, 1/2" (12.7 mm) Arbor

Application	Conc.	Grit size	Thickness	Catalog Number	Figure
4-inch CBN wafering blade (medium grit, high concentration)	High	Medium	0.012"	WCBN-0045	
5-inch CBN wafering blade (medium grit, high concentration)	High	Medium	0.015"	WCBN-0055	
6-inch CBN wafering blade (medium grit, high concentration)	High	Medium	0.020"	WCBN-0065	

CUTTING FLUID

Description	Quantity	Catalog Number	Figure
Water-based wafer cutting fluid	16 oz (0.47l) 32 oz (0.94l)	WL-3000-16 WL-3000-32	
Oil-based wafer cutting fluid	16 oz (0.47l) 32 oz (0.94l)	OL-3000-16 OL-3000-32	

Dressing stick.

Description	Quantity	Catalog Number	
Dressing sticks (medium grit)	each	DRES-0050	

-Chapter 4 Specimen mounting 에서 계속-